#1 brand among doctors!
最多医生选用品牌¹

1st Mom's choice
Leading brand²

New

AgeSmart™ for the whole family

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Infrared ear thermometer
博朗耳温计

IRT 6520 / IRT 6030


Do you know the definition of fever changes with age? My thermometer does.

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Clinical research shows that the definition of fever changes with a child's age³. Many parents and guardians are not aware of this. What is fever in a newborn can be a normal temperature in a four-year-old. To help them interpret a child's temperature at first sight, Braun has developed unique AgeSmart™ technology.

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Braun AgeSmart™ technology provides age-adjustable fever guidance to take the guess work out of interpreting temperature for the whole family.

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Fever guidance 发烧指示功能

Large color-coded display provides instant fever detection for the whole family.

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More accurate 测温准确

Exclusive pre-warmed tip technology is designed to minimize the cooling effect of the tip inside the ear canal, providing accurate measurements time after time.

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ExacTemp™ technology uses a light and beep to confirm proper positioning.

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Night light feature 夜光功能

Helps take a temperature even when it's dark.

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Memory function 记忆功能

-IRT 6520: Recalls last 9 temperature readings.

-IRT 6030: Recalls last temperature reading.

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Safe and hygienic 使用安全及卫生

Braun disposable lens filter to eliminate the risk of cross-contamination and ensure accurate readings.

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Clinically proven more accurate than rectal or temple⁴.

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
Remove the thermometer from the protective case.

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

Fit the probe snugly into the ear canal.

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Press and release the Start button.

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
Temperature displays.

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(IRT 6030)


Age range 年龄组别	Green Normal temperature 绿色 体温正常	Yellow Moderate fever 黄色 轻度发烧	Red High fever 红色 高烧
0-3 Mos 0-3 个月	> 35.8 – 37.4°C		> 37.4°C
>3-36 Mos 3-36 个月	> 35.4 – 37.6°C	> 37.6 – 38.5°C	> 38.5°C
>36 Mos adult 36 个月以上 - 成年人	> 35.4 – 37.7°C	> 37.7 – 39.4°C	> 39.4°C

Age-adjustable fever guidance for the whole family. (IRT 6520)

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• Always replace disposable lens filters to maintain accuracy and hygiene.

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What is fever? 何谓发烧？

If your body temperature is consistently above your normal temperature range*, you are said to have fever. Fever is a defense mechanism the body uses to fight infections. However, fever alone isn't an illness, it's a symptom.

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* Practice taking temperatures on yourself to determine the normal temperature range.

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Would the tip on the thermometer affect accuracy reading? 测温头会否影响体温的准确度？

Yes. The tip on other ear thermometers can cool the area where the reading is taken, leading to a higher risk of inaccurate readings. The patented Braun pre-warmed tip is designed to minimize the cooling effect of the tip inside the ear canal, providing accurate measurements time after time.

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How can parents reduce the risk of infection when their children have a fever? 小孩发烧时，父母有甚么方法能减低交叉感染的机会？

Infants and children can suffer from fever many times throughout the year, and parents usually take their child's temperature several times during each bout of fever, depending on how unwell the child feels.

Taking measurement with Braun disposable lens filter enable hygienic use without having to clean or disinfect the sensor and virtually eliminate cross-contamination. What's more, fitting a new Braun lens filter each time ensures that dirt or ear wax don't build up and clog the filter membrane, which could lead to inaccurate temperature readings.

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These products are manufactured by Kaz Europe Sàrl under license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

ThermoScan is a registered trademark of Kaz Europe Sàrl.

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© 2015, Kaz Europe Sàrl. All rights reserved.

Watch product video
看产品影片