

Mini Centrifuge



Model MCF-2360

Features

- Compact Design
- Clear PC Lid and ABS body resists breakage
- Priced as a personal centrifuge
- Flexible as a multi-purpose microcentrifuge
- Includes all adapters and accessories for
 - 1.5/2.0 ml microcentrifuge tubes
 - 0.4/0.5 ml PCR tubes
 - 0.2 ml PCR tubes
- Additional switch for better control
- Auto-start with auto-timer
- Auto-stop when lid is open during runs
- Unique anti-slip design

Specifications

Speed (rpm/g)	6,600/2,200
Capacity	1.5/2.0 ml x 6 tubes 0.4/0.5 ml x 6 tubes 0.2 ml x 2 strips of 8 tubes
Dimension (mm)	160(W) x 161.5(D) x 102(H)
Housing Color	Black/Light Gray
Power	100-220VAC 50/60Hz
Weight	1.1kg

Mini Centrifuge with Timer



Model MCT-2370

Features

- Compact Design
- 1-10 min. Timer
- Clear PC Lid and ABS body resists breakage
- Priced as a personal centrifuge
- Flexible as a multi-purpose microcentrifuge
- Includes all adapters and accessories for
 - 1.5/2.0 ml microcentrifuge tubes
 - 0.4/0.5 ml PCR tubes
 - 0.2 ml PCR tubes
- Additional switch for better control
- Auto-start with auto-timer
- Auto-stop when lid is open during runs
- Unique anti-slip design

Specifications

Speed (rpm/g)	6,600/2,200
Capacity	1.5/2.0 ml x 6 tubes 0.4/0.5 ml x 6 tubes 0.2 ml x 2 strips of 8 tubes
Timer	1-10 min.
Dimension (mm)	160(W) x 161.5(D) x 102(H)
Housing Color	Black/Light Gray
Power	100-220VAC 50/60Hz
Weight	1.1kg