

A0639

MICROSCOPE SME-F (640x)



MICROSCOPE SME-F STANDARD SPECIFICATION

COMPONENT	SPECIFICATION SME-F
Eyepiece	W.F. 10 x with pointer W.F. 16 x
Objective	4 x 10 x 40 x (spring)
Stage	Mechanical Stage 125 x 110mm Movement-X : 60mm Movement-Y : 30mm
Condenser	ABBE NA 1.25
Illuminator	Build in 20W tungsten light bulb with condenser blue filter.
Packing	Outer Packing : Card Board Box Inner Packing : Protected Styrofoam Case
Accessory	Dust Cover



Note : For further information, please contact your nearest dealer.