

Microscope XSP-12

A0643



STANDARD SPECIFICATION

COMPONENT	SPECIFICATION XSP-12
Objective	10 x
Eyepiece	4 x NA 0.10 10 x NA 0.25 40 x NA 0.65
Stage	Plane with slide clips
Condenser	Disc diaphragm
Mirror	Plan - concave
Packing	Wooden Box

This Microscope is designed in reasonable structures and easy to be operated. They are of separate coarse and fine focussing knobs and fixed with focussing stop to protect the objects and slides against damage while tube being downward.

Note : For further information, please contact your nearest dealer.