



Digital Weight Scale

MODEL **HN-286**

Ultra Slim Digital Weight Scale



- 4 Sensor Accuracy Technology
- Automatic ON/OFF Function
- Ultra Slim Design
- Precise body weight management
- Tempered glass for safety
- Precision display to the closest 100g

Specifications:

Model	Digital Weight Scale HN-286
Weight Display	5kg - 180kg
Minimum Displaying Unit	0.1kg
Power	One 3V Lithium battery (CR2032 type)
Battery Life	About 1 Year (measured 4 times a day)
Weight	About 1.5kg (including battery)
Dimensions	(L) 300mm x (W) 300mm x (H) 18mm
Accessories	Instruction manual, a trial battery

A Good Sense of Health

 **OMRON HEALTHCARE Co.,Ltd**
24 Yamanouchi Yamanoshita-cho.
Ukyo-ku, Kyoto 615-0084 Japan

Asia Pacific HQ
OMRON HEALTHCARE SINGAPORE PTE LTD
438A Alexandra Road, #05-05/08,
Alexandra Technopark, Singapore 119967
www.omron-healthcare.com.sg

Distributed by:

